

Wireless Nurse Call System

Model: FBXWNCS10-LCD

Technical Product Specification Document



FORBIX SEMICON® presents a complete 10-Bed Wireless Nurse Call System. This advanced healthcare communication system provides instant, reliable wireless alerts between patients and the nurse station, ensuring prompt response and enhanced patient safety.

- **Complete 10-Bed Setup:** Everything required for ten patient beds in a single package.
- **Smart & Customizable:** Features an LCD display, programmable room names (e.g., "Room 101", "ICU-A"), and wireless configuration via WiFi.
- **Plug-and-Play Installation:** Fully wireless operation between beds and the central station means no complex wiring.
- **Professional & Reliable:** Built with the latest RF technology, backed by one year of warranty and free technical support.

Components included in the set

The model **FBXWNCS10-LCD** system includes the following core components:

- **10 × Wireless Patient Remote Controllers** (each with an A23 battery)
- **1 × Central Receiver Unit** (with a VGA female connector for monitor connection)
- **1 × Power Adapter** (110-240V AC to 12V DC)

Technical Data

Index (Table of Contents)

Wireless Nurse Call System.....	Page 2
Patient Call Button Remote Controller.....	Page 3
Nurse Station Receiver Console with VGA.....	Page 4
Regulated AC Power Supply Adapter.....	Page 5
Remote Control Battery.....	Page 6
Application Images of Wireless Nurse Call System.....	Page 7
Nurse Station Console Parallel Monitor Screen.....	Page 8
Receiver Console Web Interface Settings.....	Page 11
Revision History.....	Page 12
Company.....	Page 13

Wireless Nurse Call System, FBXWNCS10-LCD : Technical Product Document

Product designed and manufactured by **FORBIX SEMICON®**

Patient Call Button Remote Controller

Model: FBX521



Features:

3 buttons

RED light indicator, glows when any button is pressed

Top button (ON) => send ON/CALL signal to the receiver

Bottom button (OFF) => send OFF/CANCEL signal

Center button (optional), function of which can be altered as per application requirement

Technical details:

Transmitter PA +10dbm

Operating frequency: 315 – 480 MHz

System operates works on signal hop mechanism. Signal transmission happens at more than 12-15 predefined frequencies

Modulation type: GFSK

Radio IC: CC1100 (Chipcon/TI) sub-GHz frequency

Data packet error correction mechanism: FEC

Precautions and mounting locations:

Remotes can be hand operated or mounted on the a door / wall / bed / table or any flat surface

Easy mounting of remote jacket with 2 screws

For indoor purpose use only

Keep away from metallic enclosures

Do not connect the battery in opposite direction

Range of communication:

30-40 meters indoors (depending on environmental condition)

100-150 meters outdoors (in open line of sight)

Pressing the button for at least 2-3 seconds improves the range of communication

Mechanical:

Enclosure type: ABS plastic (color white)

Remote controller size: 120 x 40 x 20 mm

Remote controller weight: 45 gms

Environmental Requirement:

Operating Temperature: 0°C to 40°C

Storage Temperature: 5°C to 25°C

Relative Humidity: 5% (0°C) ~ 95% (40°C) RH

Battery and operating voltage:

12 Volt (DC battery)

Battery model: 23A

Battery type: Alkaline (dry)

Battery replaceable: Yes

Battery rechargeable: No

Battery mounting: Back cover

Battery life: 3-4years, depending on the usage

RF antenna:

PCB trace antenna / wire coil antenna / metal antenna (depends on environment necessity)

Mounting provisions:

Remote controller can be carried in hand or mounted on the wall using a jacket provided along with the package

Application and features:

Nurse call system for hospitals: remote with patients

Steam bath sauna rooms / spas / restaurants etc.

OPD token system for dispensaries: remote with doctors

Office peon calling system: remote with clerk / manager

Guest house attendant call unit: remote with customers in room

Nurse Station Receiver Console with VGA

Model: FBXLV01

Technical Product Specification



Features:

- Centralized station receiver with LCD monitor
- WiFi control to configure settings
- Directly connected to a VGA monitor (without computer)
- 3 themes: Preloaded / Appearing / Only Numbers
- Displays 20 numbers at a time (in 5×4 matrix)
- Configurable alpha-numeric room/bed numbers
- Inbuilt RTC (Real-Time Clock)

- Runs on 220AC/12VDC adapter power supply
- VGA output for external monitor connection
- Shows real-time call alerts and patient locations for efficient response tracking

- WiFi 802.11 b/g/n support
- Capable to extend display on another computer via WiFi Connection (over Hotspot Access Point or existing WiFi network)

Precautions and mounting locations:

Device kept on the a table or can be mounted on the top of door of hospital wards / office / patient rooms / guest house rooms / hotel suites etc.

For indoor purpose use only

Keep away from metallic enclosures

Range of audible beeping sound:

Around 10-12 meters (85 decibels buzzing sound on alert)

Mechanical:

Enclosure type: ABS plastic (color white)

Receiver unit size: 200 x 150 x 70 mm

Receiver unit weight: 220 gms

Environmental Requirement:

Operating Temperature: 0°C to 40°C

Storage Temperature: 5°C to 25°C

Relative Humidity: 5% (0°C) ~ 95% (40°C) RH

Technical details:

Receiver indicator sensitivity: -112dbm

Operating frequency: 380 – 480 MHz

System operates works on signal hop mechanism. Signal transmission happens at more than 12-15 predefined frequencies

Modulation type: GFSK

Radio IC: CC1100 (Chipcon/TI) sub-GHz frequency

Data packet error correction mechanism: FEC

Operating voltage:

110-220VAC/50-60Hz to 12VDC converter adapter (supplied)

Users can also plug to any 12V external battery. Active working hours depends on the battery power

RF antenna:

PCB trace antenna / wire coil antenna / metal antenna (depends on environment necessity)

Note:

The LCD monitor is not a part of the package

Regulated AC Power Supply Adapter

Model: FBXADP12V2A



Important Note:

This product (Model: FBXADP12V2A 12V 2A Power Supply Adapter) is sourced from a qualified vendor. FORBIX SEMICON® performs quality verification at the batch level to ensure baseline specifications are met prior to distribution. However, as a non-manufactured part, this adapter is expressly excluded from the standard FORBIX SEMICON® product warranty.

Environmental Requirement:

Operating Temperature: 0°C to 40°C

Storage Temperature: 5°C to 25°C

Relative Humidity: 5% (0°C) ~ 95% (40°C) RH

Mechanical:

Enclosure: 80mm (L) × 50mm (W) × 30mm (H) Input

Connector : 2 Pin Input plug

Out Cord: UL1185 18AWG VW-1 80° 300V, FT1 Line with word on the cord Length: 0.6m - 1.2m

DC Plug : 5.5mm × 2.1mm × 12mm

Input Characteristics:

Rated Voltage: 100 – 240V AC Variation

Range: 90 – 264V AC

Rated Frequency: 50/60Hz

Variation Frequency: 47 – 63Hz

Input Current: 0.6A Max. @ Any input AC voltage & output full load

Inrush Current: 35A Max. Cold start @ 240V AC input, rated output load & 25°C ambient

AC Leakage Current: 0.25mA Max. @ 240V AC input

Output Characteristics:

Voltage: 12V

Min. Load: 0.01A

Rated Load: 2A

Output Power: 24W

Combined Load / Line Regulation:

Voltage: 12V

Min. Load: 0.01A

Rated Load: 2A

Line Regulation: ± 1%

Load Regulation: ± 3%

Ripple and Noise :

Under normal voltage & nominal load, the ripple & noise are as follows when measure with Max. Bandwidth of 20MHz & parallel 10µF electrolysis capacitor & 0.1µF ceramic capacitor crossed connect at testing point.

Ripple & Noise (Max.): 120mVp-p

Turn on Delay Time: 2S Max. @ 115V AC input and output Max. Load

Rise Time: 40ms Max. @ 115V AC input and output Max. Load

Hold Up Time: 5msec Min. @ 115V AC input and output Max.

Load Efficiency: Average efficiency 87.4% minimum @ 25%, 50%, 75% & 100% of full-loading and 230V AC input (After warm up 30 minutes).

Remote Controller Battery

Model: A23



Battery type: Alkaline

Composition: Manganese Dioxide (mercury free) Battery model: A23

Chemical System: Zn / KOH-H₂O / MnO₂ Conditions

Nominal voltage: 12 V

Open circuit voltage: 12.3...12.9 V

New battery 12.2...12.9 V

Voltage drop 10.8...10.5 V after 1 year storage at 20°C all measurements at 22 ± 2°C ambient

Capacity rated: 50 mAh discharge under 10kΩ load 24hours/day

End Voltage (EV): 6.0V minimum: 47 mAh discharge under 10kΩ load; 24hours/day

End Voltage (EV): 6.0V Standard discharge time nominal: 55 h at 10kΩ load new battery

EV: 6.0V MAD: (minimum average 50 h at 10kΩ load after 12 months storage at 20°C duration)

EV: 6.0V 52 h at 10kΩ load for IQC measurement EV: 6.0V (tested within 30 days after delivery)

Discharge currents standard: 0.05 mA max.

recommended current under continuous discharge: 0.5 mA max. recommended current under pulse discharge: 15 mA

Important Note:

This product (Model: A23 12V Alkaline Battery) is sourced from a qualified vendor. FORBIX SEMICON® performs quality verification at the batch level to ensure baseline specifications are met prior to distribution. However, as a non-manufactured part, this adapter is expressly excluded from the standard FORBIX SEMICON® product warranty.

Environmental Requirement:

Operating Temperature: 0°C to 40°C

Storage Temperature : 5°C to 55°C

Relative Humidity : 5% (0°C) ~ 95% (35°C) RH

≥90 % remaining capacity after 12 months storage

Storage should not be more than Max. 3 yrs.

Mechanical:


Unit Weight: 7.9 gms

Height: 27.70 – 28.50 mm


Diameter: 9.90 – 10.6 mm

Positive Terminal Head Diameter: 5.2 – 5.8 mm

Application Images of Wireless Nurse Call System



Nurse Station Console Parallel Monitor Screen



Wireless Calling Display Console				
A01	A02	A03	A04	A05
A06	A07	A08	Room9	Room10
Room11	Room12	Room13	Room14	Room15
Room16	Room17	Room18	Room19	Room20

FORBIX SEMICON®

- **Multi-Screen Output:** The called bed number (e.g., A03) is displayed in a blue highlight box on all parallel monitors.
- **Customizable Naming:** Default bed names (Room01-Room20) can be reconfigured in the WiFi settings.
- Default beds are named as Room01..Room20 (this nomenclature can be changed from the WiFi settings page)
- The parallel monitor is a responsive html webpage and can be view on any smart phone, tablet or computer screen
- **Web Server Access:** Access the console interface via:
 192.168.4.1 (when connected to the unit's hotspot).
 http://forbix-wncs.local (when connected to an existing network).

For hotspot access, connect to the WiFi network: **FORBIX_DISPLAY**

Receiver Console: Web Interface Settings

forbix-wncs.local/settings

FORBIX SEMICON® Wireless Calling System

Browser Time: 22/12/2024 8:33:15
RTC Time: 22/12/24 11:54:32

[Update RTC](#)

Change Display Names

Room1	Room2	Room3	Room4
Room5	Room6	Room7	Room8
Room9	Room10	Room11	Room12
Room13	Room14	Room15	Room16

forbix-wncs.local/settings

Local WiFi Login/Password:

WiFi Login

Password

[Show Password](#)

Plunge Timeout (sec), 0: disabled

Plunge Timeout

Escalation Timeout (sec), 0: disabled

Escalation Timeout

forbix-wncs.local/settings

Room17	Room18	Room19	Room20
Room21	Room22	Room23	Room24
Room25	Room26	Room27	Room28
Room29	Room30	Room31	Room32
Room33	Room34	Room35	Room36
Room37	Room38	Room39	Room40

Display Layout:

(only for connected VGA monitor)

- ☐ 5x4: Appearing Room Numbers
- ☐ 5x4: Visible Room Numbers
- ☐ 10x10: Canvas (numbers 1..100)

Local WiFi Preference:

- ☐ Connect to Hotspot
- ☐ Connect to existing Network

forbix-wncs.local/settings

Escalation Timeout (sec), 0: disabled

Escalation Timeout

Continuous Beep:

- ☐ Disable continuous beeping
- ☐ Enable continuous beeping

☐ Disable Save & Restore
 (improves device life)

☐ Enable Save & Restore beeping

[Save & Restart Device](#)

[Factory Reset](#)

Receiver Console: Web Interface Settings

Web-Based Management Interface

The console operates as a web server, providing access to its configuration panel through any standard web browser.

Network Access Points:

When connected directly to the unit's Wi-Fi hotspot: <http://192.168.4.1/settings>

When connected to an existing local Wi-Fi network: <http://forbix-wncs.local/>

System Configuration & Customization

Access the settings page to configure and personalize the system.

Real-Time Clock (RTC) Sync: Synchronize or reset the system's internal clock to match internet or computer time.

Bed/Station Nomenclature: Program custom identifiers for up to 40 beds or call points.

Default naming: Room01 to Room40.

Format: Alphanumeric, with a maximum of 6 characters per identifier (e.g., Room01, ICU-08, WardA3).

Access Credentials: Set a local administrator username and passcode for secure access to the hotspot mode or for automatic connection to a saved local Wi-Fi network.

Backup & Restore: The system features persistent memory. In case of a power failure, it automatically restores the queue of unattended calls that were active prior to shutdown.

Visual Display Modes (VGA Output):

Mode	Description	Max. Call Points	Visual Feedback
1. 5x4 Pop-up	Only called numbers appear dynamically in a 5x4 grid.	40	Numbers pop up when called.
2. 5x4 Visible	All room numbers (up to 20) are statically visible.	20	Called numbers are highlighted in a blue box.
3. 10x10 Canvas	A full grid of numbers 1-100 is displayed.	100	Called numbers change the color of their corresponding box

Receiver Console: Web Interface Settings

Call Management & Alert Parameters

Fine-tune the system's behavior and alerting logic.

Plunge Timeout: Defines the duration (in seconds) after which a serviced or unattended call is automatically removed from the active display queue.

Escalation Timeout: Introduces a delay (in seconds) before a new call appears on the main display queue. This allows primary staff a brief window to respond before the call is escalated and becomes visible to a wider team (e.g., maintenance).

Audible Alert Control: A master toggle to enable or disable the continuous system beeper. When enabled, the beeper sounds continuously as long as **any** call remains in the active, unattended queue.

Revision History

Version	Changes	Time	Verified	Approved
v1.0	Initial Release Copy	Feb'2026	Yes	-

Company

Regd. Office:

FORBIX SEMICON® Technologies Private Limited

Plot No. 77, JBR Tech Park, 6th Road,
Whitefield, Bangalore
India, Pin – 560066

Legal details:

CIN:U31909KA2021PTC156093

PAN:AAECF8700F

HSN:8542 9000

GSTIN:29AAECF8700F1ZO

Trade Mark No.:4216274, Class 9

Subsidiary Division:

FORBIX OÜ

Harju maakond, Tallinn,
Kesklinna linnaosa,
Sakala tn 7-2, 10141
Estonia

Legal details:

Registry Code: 14366377

VAT No: EE102037683

Website: <https://www.forbixindia.com>

Email: sales@forbixindia.com